

# Introduction to Whole Genome Sequencing

## IVDP: Integrated Variant Discovery Pipeline



Rodrigo Pelicioni Savegnago  
Department of Animal Science

[www.github.com/rodrigopsav/ivdp](http://www.github.com/rodrigopsav/ivdp)

## IVDP: Integrated Variant Discovery Pipeline

---

### User manual and guide

---

IVDP is a collection of Bash and R scripts developed for calling variants purposes - SNPs (single-nucleotide polymorphisms) and Indels (insertions and deletions) - from Whole Genome Sequence (WGS) and RNA Sequence (RNAseq) data. It can also do gene counts of RNAseq data. IVDP combines more than 30 programs and R packages for data manipulation and bioinformatic analysis. It is optimized to run in a local server machine or HPC with slurm scheduler.

[Citation](#) [Running IVDP first time: quick steps](#) [Download IVDP](#) [Install IVDP dependencies](#)

[Using conda environments without IVDP](#)

[Running IVDP](#)

[Killing IVDP analysis](#)

[Checking main log](#)

[IVDP parameter file](#)

[Slurm parameter file](#)

[IVDP Output](#)

[IVDP Examples](#)

- Collection of Bash and R scripts;
- It integrates more than 30 programs and R packages;
- More than 50 different functions implemented;
- Install dependencies using conda:

**`./install_ivdp_dependencies.sh -f path/to/install/dependencies/programs`**


- Scalable for HPC (High-Performance Computing) that uses SLURM, Torque, HTCondor and other scheduler systems.

- To run in a local server:

```
./ivdpRun.sh -p parFile.txt
```

- To run on a HPCC with Slurm:

```
./ivdpRun.sh -p parFile.txt -c configSlurm.txt
```



```
##### INPUT DATA #####
|INPUT_DIR=./examples/pe_wgs1
|OUTPUT_DIR=./output
|OUTPUT_NAME=PEex1a
|ANALYSIS=1 # ANALYSIS=1 (for WGS) or ANALYSIS=2 (for RNAseq)
|EXTENSION_PE1=_1.fastq.gz
|EXTENSION_PE2=_2.fastq.gz
##### INPUT GENOME #####
|ASSEMBLY=RefGenome
|REFERENCE=./examples/ref_gen/reference.fa
|VARIANTS=./examples/ref_gen/variants.vcf.gz
|ANNOTATION=./examples/ref_gen/annotation.gtf
##### WORKFLOW STEPS #####
|STEP_QC_TRIM=yes
|STEP_ALIGNMENT=yes
|STEP_MDUP=yes
|STEP_BQSR=yes
|STEP_VAR_CALL=yes
|STEP_GVCF_TO_VCF=yes
|STEP_VCF_FILTER=yes
##### CHOOSE THE PROGRAMS #####
|ALIGNER_PROGRAM=bwa2
|CALLER_PROGRAM=gatk4GenomicsDBImport, bcftools
##### GENERAL PARAMETERS #####
|CALL_BY_CHROM=yes # Only for local machine
|BP_BY_CHROM=all # Only for HPC
|CHROM_SET=2
|MIN_READ_LENGTH=50
|MAX_READ_LENGTH=150
|MIN_DEPTH=3
|MAF=0.01
|MISSING=0.3
|COMBINE_VCF=partial
|THREADS=40 # Only for local machine
|BATCH=5 # Only for local machine
#####
```

**Block to define the input files and parameters**

**Block to define the reference genome / reference variants / reference annotation**

**Block to activate (yes) or deactivate (no) each step of IVDP**

**Block to choose aligner and variant callers**

**Block of general parameters: define trimming, VCF filtering, VCF combining, and parallelism parameters for variant caller**

##### INPUT DATA #####

|INPUT\_DIR=./examples/pe\_wgs1

|OUTPUT\_DIR=./output

|OUTPUT\_NAME=PEex1a

|LIST\_SAMPLES=none # OPTIONAL

|ANALYSIS=1 # ANALYSIS=1 (for WGS) or ANALYSIS=2 (for RNAseq)

|EXTENSION\_PE1=\_1.fastq.gz # PAIRED-END READS 1

|EXTENSION\_PE2=\_2.fastq.gz # PAIRED-END READS 2

|EXTENSION\_SE=\_.fastq.gz # SINGLE-END READS

## FILE NAMES

sample3\_S1\_L001\_R1.fastq.gz  
sample3\_S1\_L001\_R2.fastq.gz  
sample4\_S1\_L001\_R1.fastq.gz  
sample4\_S1\_L001\_R2.fastq.gz  
sample5\_S1\_L001\_R1.fastq.gz  
sample5\_S1\_L001\_R2.fastq.gz  
sample6\_S1\_L001\_R1.fastq.gz  
sample6\_S1\_L001\_R2.fastq.gz



## FILE NAMES

sample3\_S1\_L001\_R1.fastq.gz  
sample3\_S1\_L001\_R2.fastq.gz  
sample4\_S1\_L001\_R1.fastq.gz  
sample4\_S1\_L001\_R2.fastq.gz  
sample5\_S1\_L001\_R1.fastq.gz  
sample5\_S1\_L001\_R2.fastq.gz  
sample6\_S1\_L001\_R1.fastq.gz  
sample6\_S1\_L001\_R2.fastq.gz

## LIST\_SAMPLES FILE

sample3\_S1\_L001 ANIMAL1  
sample4\_S1\_L001 ANIMAL1  
sample5\_S1\_L001 ANIMAL2  
sample6\_S1\_L001 ANIMAL3

```
##### INPUT GENOME #####
```

```
| ASSEMBLY=RefGenome
```

```
| REFERENCE=./examples/ref_gen/reference.fa
```

```
| VARIANTS=./examples/ref_gen/variants.vcf.gz      # BQSR and rsID
```

```
| ANNOTATION=./examples/ref_gen/annotation.gtf    # Annotation,  
                                                    # Gene count
```

##### WORKFLOW STEPS #####

|STEP\_QC\_TRIM=yes

|STEP\_ALIGNMENT=yes

|STEP\_GENECOUNT=no      # Only for RNAseq: it requires ANNOTATION

|STEP\_MDUP=yes

|STEP\_BQSR=yes

|STEP\_VAR\_CALL=yes

|STEP\_GVCF\_TO\_VCF=yes

|STEP\_VCF\_FILTER=yes

##### CHOOSE THE PROGRAMS #####

|ALIGNER\_PROGRAM=bwa2, bowtie2

(bbmap, bowtie, bowtie2, bwa, bwa2, gsnap, hisat2, star)

|CALLER\_PROGRAM=gatk4GenomicsDBImport, bcftools

(bcftools, freebayes, freebayes-parallel, gatk4,

gatk4GenomicsDBImport, gatk4CombineGVCFs, platypus, varscan)

|GENECOUNT\_PROGRAM=none # Only for RNAseq

(none, htseq, featurecounts)

##### CHOOSE THE PROGRAMS #####

|FEATURE\_TYPE=exon    # Any feature from 3<sup>rd</sup> column of gtf file.

|MDUP\_PROGRAM=sambamba

(sambamba, markdupspark, or picard)

|BQSR\_PROGRAM=bqsrspark (bqsrspark or bqsr)

```
##### GENERAL PARAMETERS #####
```

```
|CALL_BY_CHROM=yes  # Only for local machine
```

```
|BP_BY_CHROM=all    # Only for HPCC
```

```
|CHROM_SET=29
```

```
|ADAPTER_TRIMMOMATIC=/home/work/rps/software/IVDP/program/03.qc_trim/adapters/TruSeq3-PE.fa:2:30:10:8:true
```

```
|MIN_READ_LENGTH=50
```

```
|MAX_READ_LENGTH=200
```

```
##### GENERAL PARAMETERS #####
```

```
| MIN_DEPTH=3
```

```
| MAF=0.01
```


```
| MISSING=0.3
```

```
| COMBINE_VCF=none      (none, partial, full)
```


```
| THREADS=40           # Only for local machine
```

```
| BATCH=5              # Only for local machine
```

COMBINE\_VCF=partial



COMBINE\_VCF=full





TRIMMING\_TIME=04:00:00  
TRIMMING\_CPU=16  
TRIMMING\_MEM=16G

**Read trim step**

ALIGNMENT\_TIME=04:00:00  
ALIGNMENT\_CPU=16  
ALIGNMENT\_MEM=32G

**Alignment step**

GENECOUNT\_TIME=04:00:00  
GENECOUNT\_CPU=8  
GENECOUNT\_MEM=8G

**Gene count step**

MDUP\_TIME=04:00:00  
MDUP\_CPU=4  
MDUP\_MEM=16G

**Mark duplicated reads step**

BQSR\_TIME=04:00:00  
BQSR\_CPU=8  
BQSR\_MEM=16G

**BQSR step**

CALLVARIANT\_TIME=04:00:00  
CALLVARIANT\_CPU=4  
CALLVARIANT\_MEM=16G

**Variant calling step (and gVCF mode)**

GVCF\_TO\_VCF\_TIME=04:00:00  
GVCF\_TO\_VCF\_CPU=4  
GVCF\_TO\_VCF\_MEM=16G

FILTERVCF\_TIME=04:00:00  
FILTERVCF\_CPU=4  
FILTERVCF\_MEM=8G

**VCF filtering step**

TRIMMING\_TIME=04:00:00  
TRIMMING\_CPU=16  
TRIMMING\_MEM=16G

BQSR\_TIME=04:00:00  
BQSR\_CPU=8  
BQSR\_MEM=16G

ALIGNMENT\_TIME=04:00:00  
ALIGNMENT\_CPU=16  
ALIGNMENT\_MEM=32G

CALLVARIANT\_TIME=04:00:00  
CALLVARIANT\_CPU=4  
CALLVARIANT\_MEM=16G

GENECOUNT\_TIME=04:00:00  
GENECOUNT\_CPU=8  
GENECOUNT\_MEM=8G

GVCF\_TO\_VCF\_TIME=04:00:00  
GVCF\_TO\_VCF\_CPU=4  
GVCF\_TO\_VCF\_MEM=16G

MDUP\_TIME=04:00:00  
MDUP\_CPU=4  
MDUP\_MEM=16G

FILTERVCF\_TIME=04:00:00  
FILTERVCF\_CPU=4  
FILTERVCF\_MEM=8G



## IVDP - Integrated Variant Discovery Pipeline

**DATE AND TIME:** 05/26/21 15:22:20  
**NAME OF ANALYSIS:** ex1a  
**TYPE OF ANALYSIS:** WGS  
**TYPE OF FILES:** Paired-end  
**NAME OF REFERENCE GENOME:** refGen  
**GENE COUNT PROGRAMS:** none  
**ALIGNMENT PROGRAMS:** bwa2  
**VARIANT CALLING PROGRAMS:** bcftools  
**MINIMUM READ LENGTH:** 50  
**MAXIMUM READ LENGTH:** 150  
**MINIMUM DEPTH FOR FILTERED LOCI:** 3  
**MINIMUM ALLELE FREQUENCY FOR FILTERED LOCI:** 0.01  
**MAXIMUM GENOTYPES MISSING RATE:** 0.3  
**NUMBER OF FASTQ FILES (SAMPLES):** 10  
**NUMBER OF INDIVIDUALS:** 10  
**LIST OF SAMPLES:** [click here](#)  
**PARAMETER FILE:** [click here](#)

DATA FILES	REPORTS	STATISTICS
<b>QUALITY CONTROL OF FASTQ FILES</b> <a href="#">QC AND TRIMMED FASTQ FILES</a>  <b>ALIGNMENT AND GENE COUNTS</b> <a href="#">ALIGNED BAM FILES</a> <a href="#">GENE COUNTS</a>  <b>VARIANTS</b> <a href="#">GVCF FILES</a> <a href="#">RAW VCF FILES</a> <a href="#">FILTERED VCF FILES</a> <a href="#">COMBINED SNP VCF FILES</a>	<b>QUALITY CONTROL OF FASTQ FILES</b> <a href="#">QC AND TRIMMED FASTQ FILES</a>  <b>ALIGNMENT AND GENE COUNTS</b> <a href="#">ALIGNED BAM FILES</a> <a href="#">GENE COUNTS</a>  <b>VARIANTS</b>  <a href="#">RAW VCF FILES</a> <a href="#">FILTERED VCF FILES</a> <a href="#">COMBINED SNP VCF FILES</a>	<b>QUALITY CONTROL OF FASTQ FILES</b> <a href="#">QC AND TRIMMED FASTQ FILES</a>  <b>ALIGNMENT</b> <a href="#">ALIGNMENT</a>  <b>VARIANTS</b>  <a href="#">RAW VCF FILES</a> <a href="#">FILTERED VCF FILES</a> <a href="#">COMBINED SNP VCF FILES</a>

Working with IVDP

[www.github.com/rodrigopsav/ivdp](https://www.github.com/rodrigopsav/ivdp)